

Overview

Planktonic ciliates are microplankton (20 - 200 μm in size; Sieburth et al. 1978). The most common methods for sampling, fixing, and staining these organisms are briefly described and/or referenced on this page. In the data sheets we include pictures obtained using these techniques.

Inverted microscopes are commonly used to quantify and identify ciliates and other microplankton in plankton samples. Unfortunately, with routine preservation methods (e.g. Lugol's iodine), identification of ciliates is often difficult. The specialised technique of silver-staining, with protargol, reveals diagnostic features of ciliates and is normally used for taxonomic descriptions. However, protargol staining is elaborate and time consuming. Therefore, it is rarely used in routine ecological studies (Montagnes & Lynn 1993).

Plankton samples preserved with Lugol's (2 to 10 %) have low loss of cell numbers, but it is important to note that this fixative shrinks cells. In the data sheets, the average biovolume provided for taxa is based on Lugol's fixed material. A conversion factor commonly used for aloricate ciliates is 190 fg C μm^{-3} , when the cells are preserved with 2 % acid Lugol's solution (Putt & Stoecker 1989). For other concentrations or fixatives, this factor may need to be adjusted (Jerome et al. 1993; Stoecker et al. 1994). Additional information on staining, fixation and carbon conversions can be obtained from manuals for microzooplankton sampling, fixation, and staining, listed in the references at the end of this page (e.g. Gifford & Caron 2000).

Plankton sampling, storage, and enumeration

Large volumes of seawater can be collected in routine, oceanographic sampling bottles (e.g. Niskin bottles). Immediately after sampling, subsamples of a few millilitres to several litres can be taken and preserved with Lugol's iodine (see below). Sampling with bottles minimises the loss of cells and cell disruption due to the turbulence and pressure caused by pumps and plankton nets (Gifford & Caron 2000).

After Lugol's fixed samples have been collected, they can be stored in a cool, dark place for weeks. NOTE - there has been little evaluation of the effects of

extended storage (Gifford and Caron 2000). Ciliate abundance can be obtained by settling the fixed samples in settling chambers and examining them under an inverted microscope (Hasle 1978, Gifford & Caron 2000).

Two sources for settling chambers used for inverted microscopy are:

<http://www.dunc-assoc.demon.co.uk>
Duncan & Associates, Jeeves Bank, Fernleigh Road, Grange-over-Sands, Cumbria LA11 7HT, UK

Email: email info@dunc-assoc.demon.co.uk, Telephone 00 44 015395 33857, Fax: 00 44 015395 34963

<http://www.kc-denmark.dk>
KC Denmark - Research Equipment, Holmbladsvej 19 - DK-8600 Silkeborg - Denmark
Tel. +45 86828347 - Fax +45 86824950, E-mail: kc@kc-denmark.dk

Lugol's fixation



Microplankton samples should be fixed immediately to avoid loss of cells. Samples should also always be rapidly added to the concentrated fixative, diluting them to the final concentration. For our samples we used a final concentration of 2 to 5 % Lugol's iodine (vol:vol).

Acid Lugol's solution (Thronsdén 1978):

Dissolve 100 g KI in 1 L of distilled water. Dissolve 50 g iodine (crystalline) in 100 ml glacial acetic acid. Mix these two solutions. Remove any precipitates. Store in the dark.

Advantages and disadvantages of Lugol's:

Lugol's acts as a fixative-preservative-stain and is better for accurately quantifying ciliates than many aldehyde-based fixatives (Stoecker et al. 1994; Thronsdén 1978). However, Lugol's masks chlorophyll fluorescence, which may be needed to recognise mixotrophic species (Gifford & Caron 2000). Furthermore, Lugol's will dissolve hard structures such as coccoliths and diatom frustules and, therefore, is not ideal for long-term storage of many plankton taxa (Gifford & Caron 2000).

Lugol's solution is also iodine-based and is relatively harmless compared to aldehyde-based or other more toxic fixatives. Furthermore, iodine enhances the sinking of cells in settling chambers. Lugol's not only fixes cells but also stains them a dark brown colour. This simplifies counting but obscures some of the characteristic features of the ciliates (e.g. macronucleus). Darkly stained specimens can be cleared with sodium thiosulphate (see DAPI staining below).

Lugol's does not necessarily preserve the cell shape and size of live specimens. Thus, comparison of live and Lugol's fixed material is not always possible. However, Lugol's fixed material can be processed in several ways: SEM (Montagnes & Taylor 1994), DAPI (see below), protargol staining (Montagnes & Lynn 1993). Thus, Lugol's is a relatively harmless and versatile fixation method, which we recommend for routine sampling of planktonic ciliates.

DAPI staining of Lugol's fixed samples



DAPI is a fluorescent dye that binds to nuclear DNA. The fluorochrome is excited by UV-light and can be observed with epifluorescence microscopy, using an inverted microscope. In conjunction with E. Lessard (Oceanography, Univ. Washington) we are developing a technique to stain Lugol's fixed plankton samples with DAPI. This allows the nuclear shape of microplankton to be used as a diagnostic feature in routine plankton analysis. This method is being prepared for publication (Strüder-Kypke et al. in prep); we provide a preliminary protocol below.

DAPI protocol:

1. Prepare a DAPI (4',6-diamino-2-phenylindole dihydrochloride) stock solution of $150 \mu\text{g ml}^{-1}$.
2. Dilute this stock solution with distilled water to a DAPI working solution of $15 \mu\text{g ml}^{-1}$. Both the stock and the working solution should be stored in the dark at $\sim 4^\circ\text{C}$.
3. Settle 2 % Lugol's fixed material in a 10 ml settling chamber.
4. With a pipette, add a small amount of saturated sodium thiosulphate (150-200 μl) to the sample. Sodium thiosulphate bleaches the sample, reducing the darkness of the Lugol's fixed cells and enhancing the brightness of the fluorescence. The quantity of sodium thiosulphate can be varied to achieve the desired effect. NOTE - using large amounts of sodium thiosulphate can bleach the organisms too much; they become colourless and difficult to find.
5. With a micropipette, add 150 μl of the DAPI working solution to the sample (final concentration: $2.25 \mu\text{g } 10 \text{ ml}^{-1}$). At this concentration, DAPI should stain the nuclei. If higher concentrations are used, DAPI may stain the entire cytoplasm.
6. Incubate the sample for 3 - 5 min in the dark.
7. The sample can be observed with epifluorescence microscopy, using an UV filter (excitation 365 nm) The stained nucleus will fluoresce bright blue. NOTE - Samples stained with DAPI should be kept away from bright light, as DAPI is light sensitive and fades quickly.

Protargol staining of Lugol's fixed samples




Silver staining with protargol (silver-proteinate) reveals diagnostic features of ciliates (e.g. the nuclei and infraciliature of both the cell body and oral regions). Through a series of steps, the silver-proteinate is preferentially bound to these structures. Then, the silver is "developed", much like the process of black & white photography, and the structures are revealed. Protargol is routinely used for taxonomic descriptions, but can also be applied to ecological studies if a detailed knowledge of taxa is required. For more information about protargol staining, see Montagnes & Lynn (1993, 1987) or the review from Foissner (1991).

Two protocols are recommended for protargol staining of Lugol's fixed samples. The first protocol, after Wilbert (1975), involves manipulation of single cells. The cells are stained in embryological dishes and the liquids are removed with micropipettes. This requires some practice but provides good stains, as it allows you to watch the developing process. The second method is quantitative (QPS - Quantitative Protargol Staining: Montagnes & Lynn 1993, 1987). In this method, the cells are concentrated on filters and embedded in agar. All steps are then performed in small staining jars. Using QPS, the ciliates of a subsample can be identified to the species level and counted; thus species abundance can be determined. Unlike Lugol's, protargol

has the added benefit of providing a permanent record.

NOTE - Lugol's fixed samples must be postfixed with concentrated Bouin's (5 % final concentration) before they can be processed (Montagnes & Lynn 1993).

SEM



Text to follow

References

Foissner W (1991) Basic light and scanning electron microscope methods for taxonomic studies of ciliated protozoa. *Europ J Protistol* 27: 313-330

Gifford DJ, Caron DA (2000) Sampling, preservation, enumeration and biomass of marine protozooplankton. In: Harris RP et al. (eds) ICES Zooplankton methodology manual. Academic Press, London, p 193-221

Hasle GR (1978) The inverted microscope method. In: Sournia A (ed) *Phytoplankton manual*. UNESCO Paris, p 88-96

Harris RP, Wiebe PH, Lenz J, Skjoldal HR, Huntley M (2000) ICES Zooplankton methodology manual. Academic Press, London

Jerome CA, Montagnes DJS, Taylor FJR (1993) The effect of quantitative protargol stain and Lugol's and Bouin's fixatives on cell size: a more accurate estimate of ciliate species biomass. *J Euk Microbiol* 40: 254-259

Kemp PF, Sherr BF, Sherr EB, Cole JJ (1993) *Handbook of methods in aquatic microbial ecology*. Lewis Publishers, Boca Raton

Montagnes DJS, Lynn DH (1987) A quantitative protargol stain (QPS) for ciliates: method description and test of its quantitative nature. *Mar Microb Food Webs* 2: 83-93

Montagnes DJS, Lynn DH (1993) A quantitative protargol stain (QPS) for ciliates and other protists. In: Kemp PF et al. (eds) *Handbook of methods in aquatic microbial ecology*. Lewis Publishers, Boca Raton, p 229-240

Montagnes DJS, Taylor FJR (1994) The salient features of five marine ciliates in the class Spirotrichea (Oligotrichia), with notes on their culturing and behaviour. *J Euk Microbiol* 41: 569-586

Putt M, Stoecker DK (1989) An experimentally determined carbon : volume ratio for marine 'oligotrichous' ciliates from estuarine and coastal waters. *Limnol Oceanogr* 34: 1097-1103

Sieburth JMcN, Smetacek V, Lenz J (1978) Pelagic ecosystem structure: heterotrophic compartments of the plankton and their relationship to plankton size fractions. *Limnol Oceanogr* 23: 1256-1263

Stoecker DK, Gifford DJ, Putt M (1994) Preservation of marine planktonic ciliates: losses and cell shrinkage during fixation. *Mar Ecol Prog Ser* 110: 293-299

Thronsen J (1978) Preservation and storage. In: Sournia A (ed) *Phytoplankton manual*. UNESCO Paris, p 69-74

Wilbert N (1975) Eine verbesserte Technik der Protargolimprägung für Ciliaten. *Mikrokosmos* 64: 171-179